

This document provides an overview of the new features that have been developed for the new easy go version 3.0 which is available from now onwards.

Most important add-on's at a glance:

1. Possibility to change Administrator password
2. No fix assignment of a USB to a Q-tag CLm doc device needed anymore
3. Simplified settings dialog
4. File Server installation possible (settings are stored only in one protected file)
5. Former settings of easy go 2.0 will be automatically converted into the new file setup of easy go 3.0



>> [Download of the Q-tag® easy go 3.0 software and user manual](#)



For questions please contact Berlinger customer support by e-mail at support@berlinger.com

Contact:

Berlinger & Co. AG
Mitteldorfstrasse 2
9608 Ganterschwil
Switzerland
www.berlinger.com

Subject: Release notes Q-tag easy go 3.0
To: easy go users
Place, date: Ganterschwil, June 19th 2017
By: Gerhard Santner